

ENODEST - INSTALLATION GUIDE

Information to consider before installing your RAYPA equipment.

INDEX

Electrical connection	Page 2
Connections graph	Page 2
Cooling water supply	Page 3
Steam generator water supply	Page 3
Drainage connections	Page 3
Included components	Page 4
Dimensions	Page 6
Environmental conditions	Page 6




OENOLOGIC DISTILLER ENODEST


ELECTRICAL CONNECTION

The following table shows the plug configuration according to international IEC and SCHUKO standards for most European Union and LATAM countries. For customers requiring other plugs and other electrical configurations, please contact our technical team at raypa@raypa.com.

MODELS	FREQUENCY	POWER	VOLTAGE	CONNECTION
ENODEST-TS	50/60 Hz	2000 W	230 (1P+N+E) V	16 A 1
ENODEST-TS-115V	50/60 Hz	2000 W	120 (1P+N+E) V	32 A 4



CONNECTIONS GRAPH



COOLING WATER SUPPLY

Decalcified water is required for cooling the equipment. Connect the COOLING WATER INLET (2) with the hose supplied* to a decalcified water supply (pressure between 1,5 and 8bar).

The water used for cooling should have a temperature equal to or lower than 25°C.

*See section of included components for more information on the technical characteristics of this hose.

STEAM GENERATOR WATER SUPPLY

Distilled water is required to operate the steam generator and is added automatically by connecting the supplied hose* to the STEAM GENERATOR WATER INLET (5) and the other end to the quick connection of the 10L tank.

*See section of included components for more information on the technical characteristics of this hose and tank assembly.

DRAINAGE CONNECTIONS

Connect the COOLING WATER OUTLET (4) of the unit with the supplied hose*. Lead the other end to a drain.

Connect the STEAM GENERATOR DRAIN OUTLET (3) of the unit with the supplied hose*. Lead the other end to a drain.

*See section of included components for more information on the technical characteristics of this hose.

INCLUDED COMPONENTS



1 reinforced NBR hose 2m long with a 3/4" threaded connection both at the equipment and a tap (gaskets included).

For:
2. COOLING WATER INLET



1 transparent silicone hose of $\varnothing 5 \times \varnothing 8$ mm and 1m long with *press-fit* connection to connect to the equipment and, at the other end, with a quick connection to connect to the supplied tank.

For:
5. WATER INLET FOR STEAM GENERATOR



2 transparent silicone hoses of $\varnothing 8 \times \varnothing 14$ mm and 1m long with a 3/8" threaded connection on one end to connect to the equipment (gaskets included) and an unthreaded end to connect to the drainage.

For:
3. STEAM GENERATOR DRAIN OUTLET
4. COOLING WATER OUTLET



1 10 liter polyethylene tank of L x D x H: 190 x 220 x 330mm with a screw lid with quick connection to store distilled water.

For:
5. WATER INLET FOR STEAM GENERATOR

INCLUDED COMPONENTS



1 glass tube for the distillation of samples with a size of 42 x 300mm.



1 reinforced glass tube for the distillation of samples with a size of 80 x 300mm.



3 glass volumetric flasks for the collection of distillates:

- 100mL of Ø x H: 61 x 170mm
- 200mL of Ø x H: 75 x 200mm
- 250mL of Ø x H: 80 x 220mm



1 anti-drip plastic tray of L x D x H: 375 x 130 x 20mm.





OENOLOGIC DISTILLER ENODEST



DIMENSIONS TO CONSIDER FOR THE INSTALLATION OF YOUR EQUIPMENT

The equipment shall be placed on a stable and flat surface suitable for the weight of the equipment. At a distance of less than 1500mm, a water supply, a drain and an electrical outlet must be provided. For safety reasons, the distance between both sides of the equipment and the wall or any other object must be 500mm and between the equipment and the rear wall it must be at least 200mm. Do not place containers, chemicals or other devices behind the equipment.

MODEL	L LENGTH	D DEPTH	H HEIGHT	HD DRAINAGES HEIGHT
ENODEST-TS	520 mm	360 mm	910 mm	100 mm
GF-10L (tank)	190 mm	220 mm	330 mm	-



OENOLOGIC DISTILLER
ENODEST

+ info



CLICK!
ACCESS THE
ENODEST
VIDEO

Find out more about our **ENODEST** distiller on our Youtube Channel

